


DrDAQ[®]

Versatile data acquisition



Data logger, oscilloscope, signal generator and more!

Use the built-in sensors for light, sound, and temperature

Measure pH using standard electrodes

Add external sensors and digital devices

Capture data from up to 20 DrDAQs on one PC

USB connected and powered

Free to download PicoLog 6 and PicoScope 6 software

Free software development kit

Example programs available to download

Free technical support

Free software updates

Compatible with Windows, macOS and Linux

Sensors and indicators

You can start using your DrDAQ Data Logger straight out of the box; it features built-in sensors for light, sound and temperature, and an RGB LED that you can program to show any one of 16.7 million colours. External sensor sockets also allow you to extend the functionality of your DrDAQ. With external sensors your DrDAQ can measure humidity, oxygen levels, external temperatures and much more. Pico Technology provide you with everything you need to know to connect, use, and even design your own sensors.


More than just a data logger

Thanks to the versatility of your DrDAQ you can also use it as an oscilloscope and spectrum analyzer. Just download and run PicoScope software from www.picotech.com/downloads and your DrDAQ becomes a single-channel scope with a 100 kHz bandwidth, 8-bit resolution and the ability to measure voltages of up to ± 10 volts.

Digital Input/Outputs

The DrDAQ includes 4 digital input/outputs. As inputs these give you even more monitoring options, and when used as outputs they enable you to use your DrDAQ to control external devices. Two of the digital I/Os include a pulse-counting function when used as inputs, and a pulse-width modulation (PWM) output capability.

But that's not all. Your DrDAQ also includes a signal generator which can be used as both a standard function generator and an arbitrary waveform generator (AWG). The AWG function enables you to create your own waveforms.




Hobbyist, student or professional: DrDAQ is all you need

The DrDAQ Data Logger has something for everyone: whether you're a teacher looking for a fun way to perform data logging experiments in the classroom, a student wanting an inexpensive introduction to data logging and oscilloscopes, a programmer who wants to use C++ to monitor and control the inputs and outputs of real-world devices, a hobbyist who wants to monitor and control their environment, or a professional wanting to measure pH under laboratory conditions – DrDAQ has everything you need.

HOBBYIST



PROFESSIONAL



STUDENT




Data logging and oscilloscope software

A PC-based data logger is only as good as its software — that's why the DrDAQ comes with a complete software package that enables you to use your DrDAQ not only as a data logger but also as an oscilloscope and arbitrary waveform generator. Pico Technology's PC-based programs are easy and intuitive to use and offer superior viewing quality, making them ideal for education and training. Data can easily be saved, printed, emailed, and exported to spreadsheets and databases.

Both PicoLog and PicoScope benefit from free upgrades, technical support, and international language versions, and are compatible with Windows, macOS and Linux.


PicoLog data logging software

PicoLog is a complete data acquisition software package for the Pico Technology DrDAQ data loggers. It provides a visual, easy-to-use interface for users to quickly set up simple or complex acquisitions and record, view and analyze data.



PicoScope oscilloscope software

PicoScope is your complete test and measurement lab in a single application. Using the PicoScope software you can use your DrDAQ Data Logger as an oscilloscope or spectrum analyzer, and control the integrated signal generator, RGB LED, and digital outputs.




PicoSDK®

Pico's software development kit, PicoSDK, is available free of charge and allows you to write your own software and interface to third-party software packages.

Pico also maintains repositories of example code on GitHub (github.com/picotech), showing how to use PicoSDK with software packages such as Microsoft Excel, National Instruments LabVIEW and MathWorks MATLAB, or with programming languages including C, C++, C# and Visual Basic .NET.

PicoSDK and the *DrDAQ Programmer's Guide* are available to download from www.picotech.com/downloads.



PicoLog software

PicoLog is a complete data acquisition software package for the DrDAQ data logger, and is fully compatible with Windows, macOS and Linux. With its clear and user-friendly layout, ideal for use with a mouse or a touchscreen, PicoLog allows you to set up the logger and start recording with just a few clicks of the mouse, whatever your level of data logging experience. Set up simple or advanced acquisitions quickly, and record, view and analyze your data with ease.

Device settings, Graph and Table

Easily set up and adjust acquisition and math channels on one or more data loggers and check their status at a glance. You can also select **Graph** view to see live data trend lines and **Table** view to see data tabulate in real-time.

Graph view

Display your data in real time, as it is collected, on up to four independent Y axes simultaneously: set them up by dragging and dropping the entries in the **Channels & Axes** panel on the right.

Give instant feedback

We want to hear from you! Click here to contact Pico with your comments.

Capture controls

Separate **Record**, **Pause** and **Reset** buttons make it harder to press any of them by mistake.

Save and Export options


Copy your graph to the clipboard, save it as a PDF, export the raw data to a CSV file, or save the data and configuration as a robust .picolog database file.

Alarms

Set up alarms to alert you to a range of events. **Alarms** can take the form of sounds, visual notifications, graph annotations and more.

Notes & annotations

Add notes about the dataset as a whole or annotations about particular points on the graph.



Pullout information panel

Manage your channel and axis settings, alarms, notes and capture information in this easy-to-read layout. Close the panel to make more room for the capture graph, and reopen it at any time.

Multiple devices

Log data on up to 20 devices at the same time. Here, two separate data loggers are in use: one DrDAQ and one TC-08 temperature data logger.

Data view

Display all the data collected so far (**Show all data**) or keep the graph scale the same and pan along as new samples appear (**Follow new data**).

Pan and zoom controls

Zoom in, zoom out, zoom to a selection or pan through the data with these tools. If you make a mistake, just click **Undo zoom**.

Cursors and annotations

Use **Cursor** to highlight the data value and time at any point on the graph, then click **Add annotation** to mark that point with a text note.

PicoScope software

The display can be as simple or as advanced as you need. Begin with a single view of one channel, and then expand the display to include any number of live channels, math channels and reference waveforms.

Auto setup button: Configures the collection time and voltage range for clear display of signals.

Channel options: Filtering, offset, resolution enhancement, custom probes and more.

Movable axes: The vertical axes can be dragged up and down. This feature is particularly useful when one waveform is obscuring another. There's also an **Auto Arrange Axes** command.

Zoom overview: Click and drag for quick navigation in zoomed views.

Trigger toolbar: Quick access to main controls.

Tools: Including serial decoding, reference channels, macro recorder, alarms, mask limit testing and math channels.


Waveform replay tools: PicoScope automatically records up to 10000 of the most recent waveforms. You can quickly scan through to look for intermittent events, or use the **Buffer Navigator** to search visually.

Zoom and pan tools: PicoScope allows a zoom factor of several million.

Signal generator: Generates standard signals or arbitrary waveforms. Includes frequency sweep mode.

RGB LED: Allows you to manually control the built-in RGB LED to show any one of 16.7 million colors.

Digital outputs: Allows you to set the characteristics of the four digital outputs on the screw terminal block of your DrDAQ.



Ruler legend: Absolute and differential ruler measurements are listed here.

Rulers: Each axis has two rulers that can be dragged across the screen to make quick measurements of amplitude, time and frequency.

Views: PicoScope is carefully designed to make the best use of the display area. You can add new scope, spectrum and XY views with automatic or custom layouts.


Automatic measurements: Display calculated measurements for troubleshooting and analysis. You can add as many measurements as you need on each view. Each measurement includes statistical parameters showing its variability.

Trigger marker: Drag the yellow diamond to adjust trigger level and pre-trigger time.

Math channels

Sometimes you need to use data from one or more measurement channels to graph and record a calculated parameter. You can use the PicoLog equation editor to set up simple math channels such as $A-B$ or more complex functions such as log, sqrt, abs, round, min, max, mean and median.

PicoLog treats math channels like any other channel, so you can still set alarms and annotate them.




Settings

Axis 1	
<input checked="" type="checkbox"/> Light CU985/190 Light	4.73
Axis 2	
<input checked="" type="checkbox"/> Sound waveform CU985/190 Sound	-0.25
Axis 3	
<input checked="" type="checkbox"/> Sound level CU985/190 dBA	57.00 dBA
<input checked="" type="checkbox"/> 50% attenuation Maths Channel	28.50 dBA

Alarms

In PicoLog, you can set up alarms to alert you to various events. These can be as simple or as complex as you like: alarms can trigger on a signal threshold or disconnection of the data logger, or you can set up a logic expression of your own. Alarms can play sounds, display visual alerts, run applications or mark when the event occurred on the graph.



Settings

Axis 1	
<input checked="" type="checkbox"/> Light CU985/190 Light	4.79
Axis 2	
<input checked="" type="checkbox"/> Sound waveform CU985/190 Sound	-0.24
Axis 3	
<input checked="" type="checkbox"/> Sound level CU985/190 dBA	57.01 dBA
<input checked="" type="checkbox"/> 50% attenuation Maths Channel	28.50 dBA

Intuitive logger and channel setup

The **Devices** view lets you set up a multichannel acquisition system in a simple way, with the option to use multiple different Pico data loggers simultaneously. PicoLog shows you an image of each connected device, so you can quickly and easily enable or disable channels and set up their properties.

On the right, you can see the device setup for the acquisition of data from a DrDAQ and one TC-08 temperature data logger.



Robust file format

At the heart of PicoLog is the file system, which stores live capture data directly to a robust database, rather than to a single file that is vulnerable to corruption and data loss. If the computer is shut down and rebooted, PicoLog will only lose the data during the outage – saving resumes when you restart the software.

This file system also means that the size of the dataset you can capture you is virtually unlimited – the only restriction is the size of your computer's hard disk!

The .pilog file format is compatible across all operating systems, and there is no need to set up a file to save to before the capture is complete. You can also save mid-capture if you wish to share the data collected so far. Since anyone can download and install PicoLog for free, you can easily share saved data with co-workers, customers and suppliers for offline post-analysis.

Data can be exported as CSV. In addition, you can export a PDF containing a graph, channel configuration, capture notes, annotation notes and alarm trigger history.

PicoSDK®

Pico's software development kit, PicoSDK, is available free of charge and allows you to write your own software and interface to third-party software packages.

Pico also maintains repositories of example code on GitHub (github.com/picotech), showing how to use PicoSDK with software packages such as Microsoft Excel, National Instruments LabVIEW and MathWorks MATLAB, or with programming languages including C, C++, C# and Visual Basic .NET.

PicoSDK and the *DrDAQ Programmer's Guide* are available to download from www.picotech.com/downloads.



Try the PicoLog software today!

PicoLog's built-in demo mode allows you to try out the full functionality of the software with a choice of virtual devices and simulated live data. You also can use PicoLog to view previously saved data, even with no device connected. Visit www.picotech.com/downloads and select **PicoLog Data Loggers** to get your copy.

Specifications

Number of channels	14
Maximum sampling rates:	
Block mode (PicoScope and PicoSDK)	1 MS/s ^[1]
USB streaming (PicoScope)	1 kS/s ^[2]
USB streaming (PicoSDK)	100 kS/s ^[3]
PicoLog maximum USB streaming rate	1 kS/s
Memory depth (block mode)	16 320 samples at 1 MS/s, 1 000 000 samples at slower rates ^[3]
Memory depth (USB streaming)	
PicoScope	1 MS
Using PicoSDK	Unlimited
^[1] Single channel only. ^[2] Per active channel. ^[3] Shared between active channels.	
Oscilloscope/spectrum analyzer channel	
Bandwidth	100 kHz
Resolution	8 bits
Input impedance	1 MΩ
Input type	Single-ended, BNC(f) connector
Input coupling	DC
Input ranges (full scale)	±1.25 V, ±2.5 V, ±5 V, ±10 V
DC accuracy	±3%
Timebases	10 μs/div to 5000 s/div
Function generator/arbitrary waveform generator (AWG) output	
Channels	1
Connector	BNC
Standard waveforms	Sine, square, triangle, ramp (up/down), DC level
Function generator frequency range	DC to 20 kHz
Signal amplitude	0 to 1.5 V in 2 mV steps
Signal offset	±1.5 V in 2 mV steps
Accuracy	±10 mV
AWG frequency range	DC to 100 kHz
AWG buffer size	4096 samples
AWG update rate	2 MS/s (maximum)
AWG resolution	10 bits

Built-in sensors/outputs				
Channel		Range	Resolution	Accuracy
Sound waveform		±100	±0.2	not calibrated
Sound level		55 to 100 dB	1 dB	5 dB
Temperature		−10 to +70 °C	0.1 °C @ 25 °C	2 °C @ 25 °C
Light		0 to 100	0.1	not calibrated
RGB LED		16.7 million colours	8 bit x 3	not applicable
Analog inputs				
Type	Connector	Range	Resolution	Accuracy
pH	BNC	0 to 14 pH	0.02 pH	Depends on sensor calibration
Redox/ORP	BNC (shared with pH)	±2000 mV @ 1 TΩ	1.2 mV	Depends on sensor calibration
Resistance	Screw terminal	0 to 1 MΩ	250 Ω @ 10 k	5%
External sensors	3 x FCC68 4/4	Measure 0 to 2.5 V (Also supply power to external sensors and identify sensor type via ID resistor.)	1 mV	1%
Configurable digital IO				
Channels	4			
Connector type	Screw terminals			
Input characteristics	0 to 5 V (Low: 0 to 0.8 V, High: 2 to 5.5 V) TTL compatible			
Output characteristics	3.3 V, 2.2 kΩ output impedance			
PWM output	Available on 2 Digital IO channels			
Period	<65535 μs			
Pulse width	≥1 μs			
Pulse counting	Available on 2 Digital IO channels			
Range	Up to 65535 counts @ 500 kHz			

Software	
PicoLog, PicoScope and PicoSDK	Available from www.picotech.com/downloads
PicoSDK example code	Available from Pico's GitHub organization page, github.com/picotech
PicoLog user interface languages	English, French, Italian, German, Spanish, Korean, Japanese, Chinese (simplified), Russian
PicoScope user interface languages	Chinese (simplified), Chinese (traditional), Czech, Danish, Dutch, English, Finnish, French, German, Greek, Hungarian, Italian, Japanese, Korean, Norwegian, Polish, Portuguese, Romanian, Russian, Spanish, Swedish, Turkish
PC requirements	
PicoLog	Microsoft Windows 7, 8 or 10, 32-bit and 64-bit versions, macOS 10.9 (Mavericks) or later, 64-bit only, Linux*, 64-bit only Hardware requirements as operating system. * PicoLog for Linux is distributed as an AppImage, so you can install it without superuser permissions: see appimage.org for further information. The software has been tested on OpenSUSE and Ubuntu.
PicoScope ^[4]	PicoScope Stable: Windows 7, 8 or 10 (32 or 64-bit). PicoScope Beta: macOS 10.9 (Mavericks) or later, 64-bit only, Linux (OpenSUSE and Ubuntu). Hardware requirements as operating system.
PicoSDK ^[4]	Only available for Windows. Drivers also available for 64-bit Linux and macOS.
PC interface	USB 2.0 full speed (USB 1.1 and USB 3.1 compatible)
^[4] PicoScope R6.11 and PicoSDK 10.6.11 are the last versions compatible with Microsoft Windows XP (SP3) and Vista SP2, and they are also compatible with the Windows versions above.	
Environmental	
Operating temperature	0 to 70 °C (20 to 30 °C for stated accuracy)
Storage temperature	–20 to +80 °C
Operating humidity	5 to 80 %RH, non-condensing
Storage humidity	0 to 95 %RH, non-condensing
General	
Additional hardware (supplied)	USB 2.0 cable and user manuals
Input overvoltage protection	±30 V
Power requirements	Powered from USB port
Dimensions	77 x 70 x 23 mm (approx 3 x 2.7 x 0.9 in) (includes BNC connectors)
Weight	60 g (approx 2.1 oz)
Compliance	European EMC and LVD standards FCC Rules Part 15 Class A RoHS compliant
Warranty	5 years

External sensors and probes

As well as the built-in sensors, DrDAQ has sockets for optional external sensors. When a sensor is plugged in the software detects it and will automatically scale the readings. For example, when a temperature sensor is plugged in the readings are displayed in °C, and for a humidity sensor the readings are displayed in % RH.

DD100 Temperature Sensor

A high-accuracy general purpose temperature sensor with a 2 metre lead. Suitable for air, surface, or liquid measurements.



Range	-10 to +105 °C (14 to 221 °F)
Resolution (at 25°C)	0.1 °C (0.18 °F)
Accuracy (at 25°C)	0.3 °C (0.54 °F)

DD011 DrDAQ pH sensor

A robust epoxy-bodied pH electrode ideal for educational use. The pH sensor consists of a standard electrode capable of measuring the full 0 to 14 pH range, and comes complete with a small bottle of storage solution to help prevent it from drying out.



Size	12 x 120 mm
Operating temperature	0 to 60 °C (32 to 140 °F)
Resolution	0.02 pH

PP163 DrDAQ humidity sensor

A sensor which measures humidity using a 'non-condensing' technique. It has a short response time and plugs into the external sensor connections of the DrDAQ.



Size	72 x 45 x 28 mm
Operating range	20 to 90 %RH
Overall accuracy	Reading ± 10%
Operating temperature	0 to 60 °C (32 to +140 °F)
Resolution	0.2 %RH
Minimum response time	60 seconds with vigorous air movement
Maximum response time	60 minutes in still air

PP214 Oxygen Sensor

A sensor which can be used to measure the percentage of oxygen in a gas.



Sensor type	Galvanic Cell (lead-oxygen with weak acid electrolyte)
Input range	0 to 100% oxygen
Accuracy (calibrated)	±3.0% over operating conditions range
Response times	< 15 seconds for 90% response < 25 seconds for 97% response
Humidity	0 to 95% Non-condensing
Operating temperature	5 to 40 °C (41 to 104 °F)
Storage temperature	-15 to 50 °C (5 to 122 °F)

TA375 Passive oscilloscope probe

Our ergonomically designed passive oscilloscope probe with switchable 1:1 or 10:1 attenuation is suitable for use with any oscilloscope with a 1 MΩ input and 10 to 35 pF capacitance range as well as the PicoScope range of USB Oscilloscopes. Passive probes don't require a power supply or batteries so are lightweight and easily portable.



	x1	x10
Attenuation	1:1	10:1
Bandwidth	DC to 10 MHz	DC to 100 MHz
Rise time	35 ns	3.5 ns
Input resistance	1 MΩ	10 MΩ ± 2%
Input capacitance	57 pF + C _S **	15 pF
Cable length	1.2 m nominal	

** scope capacitance

PP066 DrDAQ reed switch

The reed switch sensor can be used to detect the presence of a magnetic field such as from a bar magnet or an electromagnet. Alternatively, a simple external switch can be wired to the internal screw terminals. It has a fast response time of 2 ms so can be used as an alternative to a light gate for timing applications. Other uses include monitoring the amount of time a door is left open or a machine running.



PP216 magnetic induction kit

The Magnetic Induction Kit can be used to demonstrate electromagnetic induction, the direction of induced currents and variation of the size of induced EMF with velocity. It can also be used as a student-based investigation, either with full instructions or as an open investigation where the students research the concepts.



TA137 DrDAQ redox sensor

The DrDAQ redox sensor (also known as an ORP sensor) can measure redox potential in the range of -1500 mV to $+1500$ mV. Positive readings indicate an increase in oxidation state (loss of electrons), while negative readings indicate a reduction in the oxidation state (gain of electrons).



Size	12 x 160 mm
Measuring range	±1500 mV
Junction	Ceramic
Reference cell	Ag/AgCl, KCl gel
Sensitive component	Dia.6 x 2.5 mm platinum ring
Plug	BNC with 1 meter cable

Available kits and contents

The DrDAQ is available as a single unit or as part of a kit, as detailed below:

PP706 DrDAQ

- DrDAQ
- USB installation guide
- Cable: USB 2.0 1.8 m blue

PP716 DrDAQ pH logging kit

- DrDAQ
- pH electrode
- DrDAQ 2 m temperature sensor
- USB installation guide
- Cable: USB 2.0 1.8 m blue

PP707 DrDAQ Data logging kit

- DrDAQ
- pH electrode
- 2 DrDAQ 2 m temperature sensors
- DrDAQ humidity sensor
- USB installation guide
- 3 m sensor extension cable
- Scope probe x1/x10 switchable
- Cable: USB 2.0 1.8 m blue



Ordering information

Order code	Product name	Description	USD*	EUR*	GBP*
PP706	DrDAQ	USB connected data logger, oscilloscope, signal generator and more, with built-in sensors.	129	109	95
PP716	DrDAQ pH logger kit	Includes the standard pH electrode for measuring and data logging over the full pH range.	199	169	139
PP707	DrDAQ data logger kit	Includes one scope probe, two temperature sensors and one humidity sensor for additional measurements.	329	279	225



Optional accessories

Order code	Product name	Description	USD*	EUR*	GBP*
DD100	DrDAQ temperature sensor	High-accuracy general-purpose thermistor temperature sensor with a two-meter lead.	25	21	18
DD011	DrDAQ pH sensor	Epoxy-bodied standard electrode capable of measuring the full 0 to 14 pH range.	58	49	41
PP163	DrDAQ humidity sensor	Non-condensing humidity sensor with a short response time.	81	69	56
PP214	DrDAQ oxygen sensor	Measures the percentage of oxygen in a gas.	159	139	119
TA137	DrDAQ redox sensor	Measure redox potential in the range of –1500 mV to +1500 mV.	179	149	129
PP066	DrDAQ reed switch	Detect the presence of a magnetic field such as from a bar magnet or an electromagnet.	33	28	22
EL032	DrDAQ sensor extension cable, 3 m	Replacement 3 meter (about 10 ft) sensor cable for DrDAQ sensors.	17	14	11
PP216	Magnetic induction kit	The Kit provides a simple and effective method for students to investigate the concepts of Faraday's laws.	58	49	21
TA375	Passive oscilloscope probe	Passive oscilloscope probe: 100 MHz bandwidth 1:1/10:1 switchable, BNC	25	21	17
DD010	"Datalogging in Practice" by Roger Frost	This 142 page book contains a huge number of experiment ideas for teachers.	30	26	21
MI106	USB 2.0 cable, 1.8 m**	Replacement Pico blue USB 2.0 cable, 1.8 m	9	7	6
TA268	USB 2.0 cable, 0.5 m**	Pico blue USB 2.0 cable, 0.5 m	9	7	6

* Prices correct at the time of publication. Sales taxes not included. Please check www.picotech.com for the latest prices before ordering.

** Pico blue USB cables are designed and built specifically for use with Pico Technology oscilloscopes and data loggers in order to minimize voltage drop and noise. Take care to use your DrDAQ with Pico blue USB cables only.

DrDAQ science experiments

Our website has a number of suggested scientific experiments using the DrDAQ. The example below uses the addition of a pH probe to test the pH level of a range of common drinks.

Example experiment

Dentists always warn about the dangers of drinking acidic drinks. Teenagers, they suggest, who drink large amounts suffer from tooth erosion whereby the enamel covering of the tooth is dissolved by the acid in the drink. Athletes are also at risk from drinking sports drinks.

The recommended 'safe' pH level for drinks is 5.5, and anything below this can promote tooth erosion. This experiment looks at the relative pH values of different drinks and can be used to determine whether a particular drink is 'safe' according to the above criteria. It can also be used as a good introduction to further studies into tooth decay and enamel stripping.

This experiment is suitable for ages 14+ and requires some knowledge of the pH measurement.

For further information on this experiment, please visit: picotech.com/library/experiment/ph-level-drinks-drdaq.

For the full list of experiments, please visit: picotech.com/library/experiments.



UK global headquarters:

Pico Technology
James House
Colmworth Business Park
St. Neots
Cambridgeshire
PE19 8YP
United Kingdom

+44 (0) 1480 396 395
sales@picotech.com

North America regional office:

Pico Technology
320 N Glenwood Blvd
Tyler
Texas 75702
United States

+1 800 591 2796
sales@picotech.com

Asia-Pacific regional office:

Pico Technology
Room 2252, 22/F, Centro
568 Hengfeng Road
Zhabei District
Shanghai 200070
PR China

+86 21 2226-5152
pico.china@picotech.com

Errors and omissions excepted. *Pico Technology*, *PicoLog*, *PicoScope*, *PicoSDK* and *DrDAQ* are internationally registered trademarks of Pico Technology Ltd.

LabVIEW is a trademark of National Instruments Corporation. *Linux* is the registered trademark of Linus Torvalds, registered in the U.S. and other countries. *macOS* is a trademark of Apple Inc., registered in the U.S. and other countries. *MATLAB* is a registered trademark of The MathWorks, Inc. *Windows* and *Excel* are registered trademarks of Microsoft Corporation in the United States and other countries.

MM004.en-16. Copyright © 2004–2019 Pico Technology Ltd. All rights reserved.

www.picotech.com



Pico Technology



@LifeAtPico



@picotechnologyltd



Pico Technology



@picotech